

# *DIGITAL* REFRACTION UNIT

**KEPLER-17**

*LED TOUCH PANEL*

 *BLUE TOOTH*

*AUTOMATED ARM*



*Two M Ophthotronics*  
*"We are Digital"*



**Automated Arm**

# DIGITAL REFRACTION UNIT

## KEPLER-17



**LED touch panel**




**Luxurious leather finish**



**Multifunction footswitch  
(controls patient chair & phoropter arm)**



**control the unit through "APP"**



### TECHNICAL SPECIFICATIONS

Chair height Min./Max.	560 / 760mm
Chair Elevation Stroke	200 mm
Instrument Table Height	875 mm
Maximum Area Required	2120 x 1700 mm in erect position
Back-rest	Motorized
Table Top Size	32" x 16"
Rotation	180 circular
Load Lifting Capacity	175 kgs.
Auxiliary Output	230v 50 Hz/ 4 connection
Voltage on table	6v/12v/15v
Power Supply	230 v / 110 v 50/60 Hz
Gross Weight	260 Kgs.

Please Note: Specification are subject to modifications, as per emerging market demands.



**Two M Ophthotronics**  
**"We are Digital"**