

OPHTHALMIC OPERATION TABLE

MO-2012 SS



Two M Ophthotronics

Hi-Tech Ophtho Equipment

OPHTHALMIC OPERATION TABLE MO-2012



TECHNICAL SPECIFICATIONS	
Height range min./ max	640 mm / 840 mm
Lifting speed	8 mm/s
Total length	1880 mm (with head)
Total width	600 mm
Base size	920 mm / 540 mm
Input Power Supply	230 V 50/60 Hz. 110 V 50/60 Hz.
Weight capacity	170 Kgs.
Gross weight	130 Kgs.

[Please Note: Specification are subject to modifications, as per emerging market demands.]

Motorised Vertical Height Adjustment ○

enables the ophthalmologist to perform from a convenient, comfortable standing or sitting position.

Optimally Positioned Head Rest ○

with manually controlled up & down motion option allows the surgeon easy access to the patient, while offering the patient comfortable head positioning. Plus, there is ample leg room, and sufficient space for equipment foot pedals.

Noiseless Dc Motor-driven ○

Adjustable Wrist Support ○

Offers surgeons convenient stability during operations.

Built-in Head Low and Head Up Position ensures safety during emergencies. ○



Two M Ophthotronics
Hi-Tech Ophtho Equipment

www.twom.co.in