



ELLEX I

Save Eye Save Life



POT - 360

OPHTHALMIC

MOBILE OPERATING TABLE

POT - 360

OPHTHALMIC

MOBILE OPERATING TABLE

- The cost-effective mobile operating table is meets all the requirements of Ophthalmologist.
- The mobile operating table POT 360 offers many extra features and a large of positioning possibilities.
- The special design of the castors enables both straight- line travel as well as free-wheeling.

FEATURES

- Easy transportation of patients
- Trendelenburg / Reverse Trendelenburg
- More leg space
- Imported DC motor driven
- Foot switch operated

TECHNICAL SPECIFICATION

HEIGHT RANGE MIN/MAX	: 680MM/880MM
LIFTED SPEED	: 8MM
TOTAL LENGTH	: 188MM/S
TOTAL WIDTH	: 600MM
BASE SIZE	: 840MM/520MM
INPUT POWER SUPPLY	: 230V 50/60Hz
	: 110V 50/60Hz
WEIGHT LIFTING CAPACITY	: 250KGS
GROSS WEIGHT	: 130KGS



Marketed By:

ELLEX I

No. 105, Ground Floor, 3rd Cross, 9th Main Road, 9th Block, Bandemaramma
Busstop Opp. Vinayaka Layout, Nagarabhavi 2nd Stage, Bangalore - 560 072.

E-mail : ellexibng12@gmail.com

